

THE MAIZER

- Small batch sheller for research locations, single ears or bulk samples
- Minimum seed damage, utilizes rubber belt for shelling
- 15 second cycle time for 10-15 ears
- 5 second cycle time for individual ears
- Excellent clean-out capabilities, no carry over between samples
- Modular design allows you to easily add and change options
- Easy to inspect for transgenic cleanout
- Easy to maintain and service
- 110v, 20 amp standard
- 220v, single phase, 50 hz, 60 hz available
- H 61", W 32", L 72"
- Weight: 750 lbs.
- Optional cob conveyor
- Optional seed transfer system





